

Sample Concentration for Volatile Organic Compounds (VOCs) Using Preparative USEPA Method 5030 and 5035 in Conjunction with Determinative Method 8260 Using Hydrogen Carrier Gas with the Agilent HydrolInert Source and Nitrogen Purge Gas with Teledyne Tekmar Atomx XYZ Purge and Trap

Tom Hartlein¹, Amy Nutter¹, Tyler Trent¹, Simone Novaes-Card²
¹Teledyne Labs, Mason, Ohio, USA; ²Agilent Technologies, Inc.

EMS 2024



Introduction

As helium supplies become scarcer and more expensive, laboratories have been seeking alternative carrier gases. Hydrogen is a low-cost, renewable gas suitable for several Gas Chromatograph/Mass Spectrometry (GC/MS) applications. The United States Environmental Protection Agency (US EPA) last revised its 8260D method, revision 4, in February 2017, which differs from its 8260A-C² methods with an updated Bromofluorobenzene (BFB) tune check acceptance criteria to help with interactions between the carrier gas and water vapor. This is a demonstration of hydrogen carrier gas for the quantitation of volatile organic compounds (VOCs) with the method requirements of US EPA Method 8260D¹ in conjunction with US EPA Method 5030³ for the aqueous samples and US EPA Method 5035⁴ for the solid waste samples. The Teledyne Labs Tekmar Atomx XYZ combined Purge and Trap (P&T) concentrator and autosampler was paired with an Agilent 7890B GC and 5977B MS and the HydrolInert ion source to demonstrate the method requirements (Figure 1).

Experimental

Working 5, 50, and 250 parts per million (ppm) or milligram per liter (mg/L) calibration standards were prepared in methanol from the following commercially available standards: 8260B MegaMix, 8260B Acetate, California Oxygenates, VOA (Ketones), 502.2 Calibration Mix, 2-Chloroethyl Vinyl Ether, and Hexachloroethane. In total, the standards contained 97 compounds.

Nine-point average response factor (%RSD) calibration curves for both water and soil were prepared from 0.5 ppb to 200 parts per billion (ppb) for all compounds. The %RSD was calculated for each compound using four internal standards. Internal and surrogate standards were prepared in methanol from commercially available standards at a concentration of 25 ppm, after which 5 microliters (µL) was then mixed with each 5 milliliter (mL) sample for a resulting concentration of 25 ppb.

Seven 0.5 ppb water standards and seven 1 ppb soil standards were prepared to calculate the MDL calculations. Seven 20 ppb water and soil standards were also prepared for the accuracy and precision calculations of the mid-point calibration check.

All calibration, MDL, and mid-point calibration check standards were analyzed with the Tekmar Atomx XYZ conditions in Table I (water) and Table II (soil). GC-MS conditions are shown in Table III.

Experimental

Table I: Tekmar Atomx XYZ Water Method Conditions			
Standby	Variable	Desorb	Variable
Valve Oven Temp	140°C	Methanol Needle Rinse	Off
Transfer Line Temp	140°C	Water Needle Rinse Vol	7.00 mL
Sample Mount Temp	90°C	Sweep Needle Time	0.25 min
Water Heater Temp	90°C	Desorb Preheat Temp	245°C
Sample Cup Temp	20°C	Desorb Temp	250°C
Soil Valve Temp	50°C	Desorb Time	2.00 min
Standby Flow	10 mL/min	Drain Flow	300 mL/min
Purge Ready Temp	40°C	GC Start Signal	Begin Desorb
Purge	Variable	Bake	Variable
Sample Equib Time	0.00 min	Methanol Glass Rinse	Off
Pre-sweep Time	0.25 min	Water Bake Rinse	1
Prime Sample Fill Vol	3.00 mL	Water Bake Rinse Vol	7.00 mL
Sample Volume	5.00 mL	Bake Rinse Sweep Time	0.25 min
Sweep Sample Time	0.25 min	Bake Rinse Sweep Flow	100 mL/min
Sweep Sample Flow	100 mL/min	Bake Rinse Drain Time	0.40 min
Spurge Vessel Heater	Off	Bake Time	2.00 min
Purge Time	11.00 min	Trap Bake Temp	270°C
Purge Flow	40 mL/min	MCS Bake Temp	180°C
Purge Temp	20°C	Bake Flow	200 mL/min
MCS Purge Temp	20°C		
Dry Purge Temp	20°C	Trap	9
Dry Purge Time	0.50 min	Chiller Tray	Off
Dry Purge Flow	100 mL/min	Purge Gas	Nitrogen

Table II: Tekmar Atomx XYZ Soil Method Conditions			
Standby	Variable	Purge	Variable
Valve Oven Temp	140°C	Purge Temp	20°C
Transfer Line Temp	140°C	MCS Purge Temp	20°C
Sample Mount Temp	90°C	Dry Purge Time	2.00 min
Water Heater Temp	90°C	Dry Purge Flow	100 mL/min
Sample Cup Temp	40°C	Dry Purge Temp	20°C
Soil Valve Temp	100°C	Desorb	Variable
Standby Flow	10 mL/min	Methanol Needle Rinse	Off
Purge Ready Temp	40°C	Water Needle Rinse Vol	7.00 mL
Purge	Variable	Sweep Needle Time	0.25 min
Pre-purge Time	0.00 min	Desorb Preheat Temp	245°C
Pre-purge Flow	0 mL/min	Desorb Time	2.00 min
Pre-heat Mix Speed	Slow	Desorb Temp	250°C
Sample Pre-heat Time	0.00 min	Drain Flow	300 mL/min
Pre-sweep Time	0.25 min	GC Start Signal	Begin Desorb
Water Volume	10.00 mL	Bake	Variable
Sweep Water Time	0.25 min	Bake Time	2.00 min
Sweep Water Flow	100 mL/min	Bake Flow	200 mL/min
Spurge Vessel Heater	100	Bake Temp	270°C
Purge Mix Speed	mL/min	MCS Bake Temp	180°C
Purge Time	Off	Trap	9
Purge Flow	Medium	Purge Gas	Nitrogen
	11.00 min		
	40 mL/min		

Table III: Agilent 7890B GC and 5977B MSD System Conditions	
Agilent 7890B GC Conditions	
Column	DB-624 Ultra Inert, 30m x 0.25 mm, 1.4µm Film, Hydrogen carrier gas - 1.0 mL/min
Oven Profile	35 °C, hold 2 min, 15°C/min to 100°C, 30°C/min to 230°C, hold 1 min, Run Time 11.67 min
Inlet	200°C, 100:1 Split, 4.715 psi
Agilent 5977B MSD Conditions	
Temp	Transfer Line 250°C; Source 250°C; Quad 200°C
Scan	Range 35 m/z to 270 m/z, Solvent Delay 0.50 min, Normal Scanning
Gain	Gain Factor 1.00, Extractor tune

Results and Discussion

Average response factor calibration curves (%RSD) for both water and soil are shown in Table IV. For US EPA Method 8260D in conjunction with US EPA Method 5030 (water) three compounds used a quadratic regression with $r^2 \geq 0.99$. All compounds for this method passed the LLOQ recalculation of the true value, with three compounds using different calibration ranges, which are denoted at the bottom of Table IV. Moreover, all compounds passed MDL and mid-point calibration check criteria. For US EPA Method 8260D in conjunction with US EPA Method 5035 (soil) two compounds used a linear regression with $r^2 \geq 0.995$. All compounds for this method passed the LLOQ recalculation of the true value, MDL, mid-point calibration check criteria. Due to space constraints, please contact the authors for method detection limit and mid-point check data.

US EPA Method 8260D BFB tune criteria passed after the system performed an autotune (Atune). Although this data was performed with the extractor tune (Etune), future research will evaluate this analysis with the Atune. US EPA Method 8260C² has an extensive list of recommended minimum relative response factor criteria for initial and continuing calibration verification; while US EPA Method 8260D displays this table, it is noted that this criterion was developed using helium carrier gas and is not appropriate for hydrogen carrier gas due to reduced response of some analytes. Future research will evaluate the minimum relative response factors with the HydrolInert source. Initial results show an increase in relative response factors in the HydrolInert source compared to the electron impact (EI) ion source, as well as, lower %RSDs of the average response factors in the calibration curve (water).

Figure 2 displays an overlay of the calibration standards, with an inset demonstrating excellent peak shape and separation across all concentrations with minimal water interference for all target compounds.




Figure 2: Overlaid TIC of calibration standards from 5-200 ppb in water and in soil

25: Ethyl tert-butyl ether 26: 2,2-Dichloropropane 27: cis-1,2-Dichloroethane 28: 2-Butanone 29: Propionitrile 30: Methyl Acrylate 31: Methacrylonitrile 32: Bromochloromethane 33: Tetrahydrofuran 34: Chloroform 35: Dibromofluoromethane (SS) 36: 1,1,1-Trichloroethane 37: Pentafluorobenzene (IS) 38: 1,1-Dichloropropene 39: Carbon Tetrachloride 40: 1,2-Dichloroethane-d4 (SS) 41: Benzene 42: 1,2-Dichloroethane 43: Isobutyl alcohol 44: Propyl Acetate 45: Ethyl Acetate 46: tert-Amyl methyl ether 47: 1,4-Difluorobenzene (IS) 48: Trichloroethylene 49: 1,2-Dichloropropane 50: Dibromomethane 51: Methyl Methacrylate 52: Isopropyl Acetate 53: Bromodichloromethane

Table IV: US EPA Method 8260D Calibration Data in Water and Soil																	
Compound	Water Calibration (0.5-200 ppb)				Soil Calibration (0.5-200 ppb)				Compound	Water Calibration (0.5-200 ppb)				Soil Calibration (0.5-200 ppb)			
	RT	Ion	RRF ±% RSD	Avg. RRF	RT	Ion	RRF ±% RSD	Avg. RRF		RT	Ion	RRF ±% RSD	Avg. RRF	RT	Ion	RRF ±% RSD	Avg. RRF
Dichlorodifluoromethane	0.828	85	9.1	0.291	0.828	85	9.4	0.060	2-Nitropropane	4.95	43	10.0	0.089	4.93	43	10.8	0.055
Chloromethane	0.929	50	10.7	0.055	0.935	50	6.7	0.140	2-Chloroethyl Vinyl Ether	5.05	63	15.6	0.194	5.04	63	9.9	0.205
Vinyl Chloride	0.984	62	9.6	0.128	0.984	62	8.0	0.187	cis-1,3-Dichloropropene	5.15	75	6.3	0.083	5.13	75	8.3	0.035
Bromomethane	1.17	94	12.6	0.194	1.17	94	7.7	0.077	4-Methyl-2-Pentanone	5.30	43	10.0	0.756	5.37	43	10.4	1.03
Chloroethane	1.23	64	8.3	0.081	1.23	64	9.0	0.616	Toluene-d8 (SS)	5.38	98	1.4	0.574	5.37	98	0.8	0.77
Trichlorofluoromethane	1.40	101	7.7	0.725	1.40	101	6.2	0.056	Toluene	5.43	91	9.2	0.233	5.42	91	11.0	0.205
Diethyl Ether	1.60	59	5.7	0.049	1.59	59	7.0	0.096	trans-1,3-Dichloropropene	5.63	75	6.7	0.177	5.61	75	6.3	0.134
1,1-Dichloroethane	1.74	61	13.2	0.046	1.73	61	6.9	0.157	Ethyl Methacrylate	5.76	69	6.9	0.239	5.74	69	5.0	0.117
Trichlorotrifluoroethane	1.76	101	10.2	0.118	1.74	101	6.1	0.082	1,1,2-Trichloroethane	5.77	97	6.2	0.974	5.76	97	7.3	0.551
Acetone ¹²	1.77	43	1.0	0.070	1.77	43	0.998	0.189	Tetrachloroethylene	5.91	166	8.7	0.231	5.90	166	9.1	0.200
Iodomethane	1.83	142	15.7	0.083	1.83	142	13.1	0.147	1,3-Dichloropropane	5.92	76	2.9	0.06	5.90	76	12.4	0.042
Carbon Disulfide	2.01	76	11.7	0.035	1.88	76	14.7	0.041	2-Hexanone	6.03	43	10.3	0.409	6.00	43	6.3	0.326
Allyl Chloride	2.01	41	11.1	0.029	1.99	41	7.3	0.028	Dibromochloromethane	6.10	129	5.4	0.107	6.09	129	3.8	0.078
Methyl Acetate	2.05	43	6.0	0.059	2.04	43	10.2	0.117	Butyl Acetate	6.18	43	6.9	0.207	6.17	43	6.1	0.164
Methylene Chloride ¹²	2.11	49	1.0	0.060	2.08	49	0.998	0.073	1,2-Dibromoethane	6.18	107	5.0		6.18	107	8.2	
Acetonitrile	2.26	41	8.3	0.017	2.23	59	10.4	0.909	Chlorobenzene-d5 (IS)	6.61	117			6.60	117		
tert-Butyl alcohol	2.26	59	6.9	0.147	2.26	41	14.9	0.609	Chlorobenzene	6.63	112			6.63	112		
Acrylonitrile	2.30	53	9.5	0.045	2.26	53	8.5	0.022	1,1,1,2-Tetrachloroethane	6.71	131	8.2	0.584	6.70	131	7.0	0.857
trans-1,2-Dichloroethane	2.34	61	9.1	0.038	2.31	61	5.3	0.084	Ethylbenzene	6.74	91	4.9	0.308	6.74	91	5.5	0.385
Methyl tert-butyl ether	2.37	73	4.1	0.670	2.34	73	4.3	0.126	m, p-Xylene	6.84	91	10.9	0.707	6.84	91	12.6	1.14
1,1-Dichloroethane	2.68	63	10.3	0.167	2.63	63	5.4	0.066	o-Xylene	7.13	91	11.2	0.549	7.13	91	13.2	0.864
Vinyl Acetate	2.75	43	13.9	0.103	2.70	43	6.3	0.09	Styrene	7.14	104	10.4	0.631	7.13	104	13.2	0.921
Chloroprene	2.76	53	12.7	0.086	2.72	53	7.9	0.128	Bromoform	7.24	173	10.8	0.685	7.24	173	11.8	0.922
Diisopropyl ether	2.80	45	5.4	0.234	2.74	45	5.6	0.223	Amyl Acetate	7.32	43	4.7	0.360	7.32	43	6.6	0.388
Ethyl tert-butyl ether	3.08	59	3.2	0.271	3.04	59	4.9	0.222	Isopropylbenzene	7.40	105	12.1	0.087	7.40	105	7.5	0.101
2,2-Dichloropropane	3.14	77	6.7	0.195	3.10	77	8.4	0.222	cis-1,4-Dichloro-2-Butene	7.44	75	11.0	1.04	7.49	75	10.5	1.43
cis-1,2-Dichloroethane	3.14	61	9.8	0.098	3.10	61	5.1	0.159	Bromofluorobenzene (SS)	7.49	95	4.9	0.421	7.49	95	3.2	0.449
2-Butanone	3.17	43	9.6	0.032	3.12	43	7.9	0.018	Bromobenzene	7.58	77	11.2	0.241	7.49	75	8.8	0.144
Propionitrile	3.18	54	11.6	0.017	3.14	54	8.5	0.011	1,1,2,2-Tetrachloroethane	7.59	83	14.7	0.126	7.58	77	10.0	0.598
Methyl Acrylate	3.26	55	9.0	0.057	3.22	55	4.6	0.049	1,2,3-Trichloropropane	7.61	75	6.3	0.150	7.59	83	6.9	0.184
Methacrylonitrile	3.31	41	13.6	0.040	3.27	41	7.5	0.026	Trans-1,4-dichloro-2-butene	7.63	75	13.1	0.069	7.63	75	9.4	0.037
Bromochloromethane	3.32	49	11.3	0.033	3.28	49	6.4	0.113	n-Propylbenzene	7.68	91	11.8	0.990	7.68	91	11.2	1.37
Tetrahydrofuran	3.38	42	10.0	0.014	3.34	42	14.6	0.011	2-Chlorotoluene	7.72	91	11.1	0.591	7.72	91	10.0	0.743
Chloroform	3.40	83	8.9	0.308	3.36	83	7.3	0.281	4-Chlorotoluene	7.79	91	10.1	0.727	7.79	91	12.3	0.946
Dibromofluoromethane (SS)	3.52	113	4.9	0.284	3.48	113	1.1	0.241	1,3,5-Trimethylbenzene	7.80	105	11.0	0.942	7.80	105	14.2	1.27
1,1,1-Trichloroethane	3.54	97	9.9	0.285	3.50	97	8.0	0.311	Pentachloroethane	8.00	117	19.5	0.055	8.00	167	6.3	0.211
Pentafluorobenzene (IS)	3.60	168			3.56	168			tert-Butylbenzene	8.01	119	9.8	1.04	8.00	119	11.6	1.45
1,1-Dichloropropene	3.67																